

# MANUAL INFORMATION SPEAKER

## WHAT'S IN THE BOX?



- Information speaker
- premium quality audio cable

## GENERAL INFORMATION

The Information Speaker has been specially developed for harmonic frequency sound therapy via the Information Medicine online platform:

**[www.informationmedicine.org](http://www.informationmedicine.org)**.

**Information medicine activates the self-healing capacity and internal intelligence in the body of both humans and animals. The harmonic tones "inform" the body with frequencies in the form of sound and light. This enables healing at the deepest level of body, mind & spirit.**

Sound therapy with harmonic frequencies is very suitable for preventive health promotion to stay fit, vital and resilient. It also offers valuable support for a variety of complaints of body, mind & spirit. For example, when recovering from physical injuries and disorders, but also for mental-emotional and energetic complaints. In our view, this relatively new method of treatment is *the* health revolution for the future.

The more natural and purer the information stimulus is, the better the body will absorb the information. The body's cells recognize the harmonics as "the original resonance of the cells" and absorb what they need to restore harmony. The quality of the speaker and of the audio cable are essential components for this form of therapy.

## DESCRIPTION INFORMATION SPEAKER

By using the Information Speaker with the supplied premium audio cable, you are assured of the optimal transfer of the information, in the form of harmonic frequencies and vibrations.

You can order a suitable premium quality audio cable from our website. We have aux cables ,et 3.5 mm jack connection , with USBC connection and with Iphone Lightning connection in stock. The audio cable provides a more stable transfer of sound information compared to wireless playback via Bluetooth or Wifi. One end of the cable plugs into the Information Speaker and the other end plugs into your device.

**The Information Speaker** is not actually a speaker, but a sound exciter designed to transmit acoustic energy (sound and vibration) through a contact surface, in this case the body. So your body actually turns into a speaker through the vibrations transmitted through the sound exciter. So the action of our harmonic frequencies is not exclusively through your ears.



# MANUAL INFORMATION SPEAKER

The Information Speaker transmits the healing information in three ways via sound frequencies, electromagnetic frequencies and through tactile vibrations.

## CONNECTING THE INFORMATION SPEAKER

Connect the Information Speaker to your playing device (PC/laptop/tablet/smartphone) with the matching audio cable. To play the Information Medicine sound files, your playing device must be suitable for online streaming. Further information on playing our sound files can be found on our website: [www.informationmedicine.org](http://www.informationmedicine.org)

## PLACEMENT OF THE INFORMATION SPEAKER

*For people:*

You place the **Information Speaker** with the vibrating side on the chest, at the height of the thymus. This is about 2 cm below the dimple of the collarbones in your neck. You can also place the Information locally in case of complaints (for example on the knee in case of knee complaints). Preferably with the vibrating plate on bare skin. Please note that the Information Speaker is not drop or impact resistant.

*For animals:*

You place the **Information Speaker** with the vibrating side at the level of the transition from the lumbar vertebrae to the sacrum (a small dimple can often be felt here in the lower back) or locally in case of complaints (for example, in the case of a tendon injury locally on the affected tendon). Please make sure that you fasten/place the Information Speaker safely and securely, as the Information Speaker is not drop or impact resistant. For horses, for example, you can use a thin cotton summer rug/fly rug and sew a pocket (for example a cotton washcloth) onto it in which you can put the Information Speaker. For very small animals, we do not recommend using the Information Speaker on the animal for practical reasons. In that case, you place the speaker in the room with the animal.

## VOLUME

It is not necessary or better to turn the volume very loud. Make sure that the harmonic frequencies can be heard and felt without being annoying or too loud. We recommend setting the volume between 25 and 65 decibels. Limit the volume to a maximum of 50-55 decibels in young children. If desired, you can use a 'decibel meter' app on your smartphone to check whether you have set the volume correctly. In that case, hold the decibel meter next to the Information Speaker during the measurement and check the volume for both the high and low tones. But most importantly, you will experience the vibration function of the sound exciter.

NOTE You cannot compare the volume and sound of the Information Speaker one-to-one with a regular speaker, which will usually sound much louder. The Information Speaker transmits the sound vibrations through the vibration plate, thus providing a tactile (tactile) sound experience.

*Volume adjustment:* You control the volume of the Information Speaker via your playing device (PC/laptop/tablet/smartphone).



# MANUAL INFORMATION SPEAKER

*Please note:* some smartphones automatically limit the playback volume, because they consider the speaker as a headphone. You can adjust this manually in the settings on your smartphone if desired.

## TECHNICAL SPECIFICATIONS INFORMATION SPEAKER

With built-in full range exciter speaker, with neodymium motor and 2.54cm voice coil and 10 watts RMS power handling

Size: Ø 6.7 cm/ 2,6 inch

Weight: 100 grams/ 3,5 oz

Total output power: 10 watts

Impedance: 4 ohms

Resonance frequency: 306 Hz

Connection: 3.5mm jack

Sound reproduction: Mono

Internal Antenna: Yes - 2.4Ghz Wifi Internal Antenna 5dBi

Amplifier built-in: No

## WARRANTY CONDITIONS

Information Medicine BV guarantees that the product is free from manufacturing and/or material defects, provided that:

- The product was sold directly to you through Information Medicine BV or by an authorized reseller.
- Was delivered to you in new condition.
- Was delivered in the original packaging.
- That the product has been used under normal circumstances.

The invoice you receive from us after placing and paying for your order is also your proof of warranty. So keep it well. If you have lost the invoice, you can find in your personal dashboard under the heading invoices.

If you do not have an account with us, please contact our customer service at [contact@informationmedicine.org](mailto:contact@informationmedicine.org).

Within the warranty period, the Information Speaker will be repaired free of charge or replaced for the same or an equivalent product if repair is not possible. No shipping costs are charged when warranty applies.

If an order is delivered visibly damaged by the relevant delivery service, refuse to accept it and do not sign for receipt. The product will automatically be returned to us and we will arrange for the delivery of a new one.

If you experience problems with an Information Speaker that you purchased second-hand, defects are not covered by the warranty. However, you are always welcome to contact us, so we can see if we can still provide a solution.



# MANUAL INFORMATION SPEAKER

**There are a number of situations that are not covered by our warranty, namely:**

- Damage or defects due to falling, impact or made intentionally.
- Defects due to improper use, for example, connecting the speaker incorrectly or damage due to overloading the speaker, etc.
- Normal wear and tear from use.
- Defects due to unauthorized third-party modifications and repairs.
- Moisture, water and fire damage.
- Damage or malfunction due to lightning strikes and electromagnetic induction due to lightning.

If you have problems with an Information Speaker that you bought second-hand, defects are not covered by the warranty, but in that case you can always contact us at [contact@informationmedicine.org](mailto:contact@informationmedicine.org) so that we can see if we can still provide a solution.

## REPAIRS

It is of course very inconvenient if a product you have purchased from us becomes defective or does not function properly (anymore). In that case, please contact our customer service at [contact@informationmedicine.org](mailto:contact@informationmedicine.org). Our customer service will advise you on the procedure to follow. If the product needs to be shipped, please do so preferably in the original packaging. Always include a copy of the invoice (your proof of purchase) and an accurate description of the complaint.

## WARRANTY PERIOD

When purchasing a product, you always have a legal guarantee. If you buy a product, you may assume that it is a sound product and that you can use it for the purpose for which it is intended. If a product does not function properly, you are entitled to free repair, a new product, or a refund. What you can expect from a product depends on the type of product, which is why there is no general legal warranty period in the Netherlands. However, within European law there is a fixed guarantee period.

- In accordance with European law, we provide a 24-month warranty on our Information Speaker, starting on the day of purchase as stated on the invoice.
- If there is a defect in the product within 6 months after purchase, the burden of proof regarding the soundness of the product lies with us, the seller.
- If there is a defect in the product after 6 months after purchase, then the burden of proof regarding the soundness of the product lies with you, the buyer/consumer.
- Note: Please note that the Information Speaker is not drop or impact resistant. Fall and impact damage is therefore not covered by the warranty. For more information about the warranty of our products, please refer to our warranty conditions. You can find it on our website: [www.informationmedicine.org/terms-and-conditions/](http://www.informationmedicine.org/terms-and-conditions/)

